## Marketed and Distributed by:





# H-IOO Urine Analyzer





# **H-100**

Urine Analyzer

#### **Product Characteristics:**

- Adopting the advanced high luminosity cold light source with 4-wavelength, which improves the sensitivity, accuracy, specificity, and reduces the interference from ambient light
- Adopting automatic waste handling system, which avoids cross-contamination between samples
- Automatically rectifies the test results influenced by non-specificity, pH, specific gravity, and color
- Built-in thermal printer with high speed and low noise; External stylus printer
- Connectable with DIRUI urine sediment analyzer
- Users can set an abnormal value flag by themselves
- International, regular and symbol system units display for option







### **Technical Specification:**

- Test items: urobilinogen, bilirubin, ketone, Creatinine, Micro-albumin, blood, protein, nitrite, leukocytes, glucose, specific gravity, pH and VC
- Test wavelength: 525nm, 572nm, 610nm, 660nm
- Test principle: Photoelectric colorimetry
- Suitable strips: DIRUI H8, H10, H11, H11-MA, H13-Cr, Urinalysis strips
- Test throughput: 120 strips/h or 60 strips/h optional
- Data storage: 1000 patient results
- Computer interface: RS-232 port; parallel printer interface
- Display: 240×64 dot-matrix LCD
- Language: Chinese, English, Russian, Polish, Italian, Spanish, Portuguese, Turkish, German, French
- Power supply: 100~240VAC, 50Hz/60Hz
- Power: 40VA
- Dimensions: 376mm×316mm×170mm
- Weight: 3.6kg
- Printer: Built-in thermal printer







GENWORKS HEALTH PRIVATE LIMITED
Unit 522-524, 5th Floor, Gamma Block, Sigma Soft Tech Park,
Bangalore - 560 066, India

Contact Us: +91 99400 84724 Ivdhelpdesk@genworkshealth.com

